



**COOLIEF* TRANSDISCAL*
COOLED RADIOFREQUENCY
Y-CONNECTOR CABLE**

Instructions for Use

Rx Only	Caution	Consult instructions for use	Dispose of properly
---------	---------	------------------------------	---------------------

HALYARD* COOLIEF* TRANSDISCAL* Cooled Radiofrequency Y-Connector Cable

Rx Only: Federal (USA) law restricts this device to sale by or on the order of a physician.

Device Description

The HALYARD* COOLIEF* TRANSDISCAL* Cooled Radiofrequency Y-Connector Cable provides a path for delivery of radiofrequency (RF) energy and temperature signals to/from the HALYARD* COOLIEF* TRANSDISCAL* Cooled Radiofrequency Probe(s). The TRANSDISCAL* Y-Connector Cable is used with one or two TRANSDISCAL* Probes, one or two HALYARD* COOLIEF* TRANSDISCAL* Cooled Radiofrequency Introducers and a HALYARD* COOLIEF* Cooled Radiofrequency (RF) Pain Management Generator or KIMBERLY-CLARK® Radiofrequency Generator (formerly Baylis Pain Management Generator), Version 2.0 or higher (PMG-115-TD/PMG-230-TD/PMG-ADVANCED).

HALYARD* COOLIEF* TRANSDISCAL* Cooled Radiofrequency Y-Connector Cable (Fig. 1)

- Shipped non-sterile and must be sterilized as per Instructions for Use prior to use.
- Used to attach the TRANSDISCAL* Probe(s) and Introducer(s) to a RF Generator.
- Equipped with three connectors:
 - 2 x 5-Pin Female RF Probe Connectors
 - 1 x 14-Pin Male HALYARD* RF Generator Connector

Indications for Use

HALYARD* COOLIEF* TRANSDISCAL* Cooled Radiofrequency Y-Connector Cables are used in conjunction with the RF Generator to create radiofrequency lesions in nervous tissue.

Contraindications

For patients with cardiac pacemakers, a variety of changes can occur during and after the treatment. In sensing mode the pacemaker may interpret the RF signal as a heartbeat and may fail to pace the heart. Contact the pacemaker company to determine if the pacemaker should be converted to a fixed-rate pacing during the radiofrequency procedure. Evaluate the patient's pacing system after the procedure.

Check the compatibility and safety of combinations of other physiological monitoring and electrical apparatus to be used on the patient in addition to the RF Generator.

If the patient has a spinal cord, deep brain, or other stimulator, contact the manufacturer to determine if the stimulator needs to be in the bipolar stimulation mode or in the OFF position.

This procedure should be reconsidered in patients with any prior neurological deficit.

The use of general anesthesia is contraindicated. To allow for patient feedback and response during the procedure, it should be performed under local anesthesia.

Systemic infection or local infection in area of the procedure.

Blood coagulation disorders or anticoagulant use.

⚠ Warnings

The TRANSDISCAL* Y-Connector Cable is shipped non-sterile and must be cleaned and sterilized as per the Instructions for Use before every procedure.

The TRANSDISCAL* Y-Connector Cable is a reusable device. Failure to properly clean and sterilize the device can cause patient infection.

The TRANSDISCAL* Y-Connector Cable must be used with the TRANSDISCAL* Probe(s).

The RF Generator is capable of delivering significant electrical power. Patient or operator injury can result from improper handling of the cable, particularly when operating the device.

Laboratory staff and patients can undergo significant x-ray exposure during radiofrequency procedures due to the continuous use of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure.

Discontinue use if inaccurate, erratic or sluggish temperature readings are observed. Use of damaged equipment may cause patient injury.

Do not modify HALYARD* Equipment. Any modifications may compromise safety and efficacy of the device.

When the RF Generator is activated, the conducted and radiated electrical fields may interfere with other electrical medical equipment.

During power delivery, the patient should not be allowed to come in contact with grounded metal surfaces.

Do not remove or withdraw the device while energy is being delivered.

⚠ Precautions

Do not attempt to use the TRANSDISCAL* Y-Connector Cable before thoroughly reading these Instructions for Use and the Instructions for Use for the TRANSDISCAL* Probe and TRANSDISCAL* Introducer and the User's Manual for the RF Generator.

Apparent low power output or failure of the equipment to function properly at normal settings may indicate: 1) faulty application of the dispersive electrode or 2) power failure to an electrical lead. Do not increase power level before checking for obvious defects or misapplication.

To prevent the risk of ignition, make sure that flammable material is not present in the room during RF power application.

Only physicians familiar with RF lesion techniques should use the TRANSDISCAL* Kit components.

It is the physician's responsibility to determine, assess and communicate to each individual patient all foreseeable risks of the RF lesion procedure.

Adverse Events

Complications with the use of this device are those associated with the procedures that utilize a RF probe.

Product Specifications

The TRANSDISCAL* Y-Connector Cable provides a path for delivery of RF energy and temperature signals to/from the TRANSDISCAL* Probe(s).

Model Number	TDX-Y-TSW-TDP
Length	8 ft.

Storage Requirements

TRANSDISCAL* Y-Connector Cable must be stored in a cool, dry environment.

Inspection Prior to Use

Perform the following checks before the patient is presented for the procedure. These tests will allow you to verify that the equipment you will use is in proper working order. Do these tests in a sterile environment.

1. Check that the TRANSDISCAL* Y-Connector Cable has been sterilized. The TRANSDISCAL* Y-Connector Cable should be sterilized before each use.
2. Visually inspect the cable for damage. Look for damage to the insulation, connectors and part labels. Ensure that there is no visible damage, such as discoloration, cracks, label fading or kinks. Do NOT use damaged equipment.
3. Ensure the TRANSDISCAL* Y-Connector Cable is dry. Residual moisture can cause improper function.

Equipment Required

RF lesion procedures should be performed in a specialized clinical setting with fluoroscopic equipment. The equipment required for the procedure is as follows:

- COOLIEF* TRANSDISCAL* Cooled Radiofrequency Introducer(s)
- COOLIEF* TRANSDISCAL* Cooled Radiofrequency Probe(s)
- COOLIEF* TRANSDISCAL* Cooled Radiofrequency Y-Connector Cable
- COOLIEF* Cooled Radiofrequency Sterile Tube Kit
- COOLIEF* Cooled Radiofrequency Peristaltic Pump Unit
- Dispersive Electrode
- Radiofrequency Generator (PMG-115-TD/PMG-230-TD/PMG-ADVANCED) Version 2.0 or higher

Instructions for Use (Fig. 2 and Fig. 3)

1. Assemble all required equipment for the intended procedure.
2. Perform a visual check on all of the equipment. Ensure all components are in good working order, and have not been damaged in any way. Do NOT use damaged equipment.

⚠ Warning: Use aseptic techniques while handling any parts of the TRANSDISCAL® Y-Connector Cable that enter the sterile field.

3. Connect the TRANSDISCAL® Y-Connector Cable to the RF Generator. When connecting the connectors, an audible click will ensure that the connection is properly secured. If connecting TRANSDISCAL® Y-Connector Cable directly to the RF Generator, connect the 14-pin male connector to its front panel. When using only the TRANSDISCAL® Y-Connector Cable, the Probe A side of the cable must be used for placement.
4. Connect the TRANSDISCAL® Y-Connector Cable to either the TRANSDISCAL® Introducer(s) or to the TRANSDISCAL® Probe(s).
5. After placement of the Introducers or after the TRANSDISCAL® Procedure, (when disconnecting the TRANSDISCAL® Introducer[s] and the TRANSDISCAL® Probe[s] from the TRANSDISCAL® Y-Connector Cable), do so by pulling on the connector body.

⚠ Caution: Prevent damage to your cable. When pulling the connectors apart, be sure to pull on the connector body and not the cable.
6. Prepare the reusable TRANSDISCAL® Y-Connector Cable for cleaning and sterilization. Transfer the used TRANSDISCAL® Y-Connector Cable to a carrying surface and cover it with a wet cloth to ensure that blood and other contaminants do not dry on the surface.

Cleaning and Sterilization Instructions

⚠ Danger

The TRANSDISCAL® Y-Connector Cable is shipped non-sterile and must be cleaned and sterilized as per these Instructions for Use prior to each use. Failure to properly clean and sterilize the device can cause patient injury and/or the communication of infectious diseases from one patient to another.

Important

The manufacturer recommends the user follow a quality control program for each sterilization cycle that meets or exceeds American Operating Room Nurses (AORN) Standards, Recommended Practices & Guidelines – 2000. This program includes, but is not limited to recording:

- Type of sterilizer and cycle used
- Lot control number
- Load contents
- Exposure time and temperature, if not provided by a recording chart
- Operator's name
- Results of sterilization process monitoring (i.e., chemical, mechanical, biological)

Cleaning and Decontamination

1. Ensure that blood and other contaminants do not dry on the TRANSDISCAL® Y-Connector Cable.
2. Rinse the TRANSDISCAL® Y-Connector Cable with deionized water until colorless run-off water occurs. Once the water runs clear soak the cable (except for the connectors) in deionized water at 22°C–48°C for 1 minute. Remove TRANSDISCAL® Y-Connector Cable from the water and scrub it with a soft bristle brush until it is visually clean.

Note: Do not let the connectors soak. Wipe connectors as necessary until they are visually clean.
3. Soak the cable (except connectors) in an enzymatic cleaning solution for 20 minutes. Ensure that the temperature of the solution is below 55°C. Scrub again with a soft bristle brush, and rinse thoroughly using deionized water until all traces of detergent residue are removed.
4. Visually inspect the cable again for debris. If any debris is present, repeat steps 2 and 3.
5. Dry the surface of the device on the outside with a clean, dry towel. Place the TRANSDISCAL® Y-Connector Cable in a Sterilization and Storage Tray.

Sterilization

The following sterilization methods have been validated for use with the TRANSDISCAL® Y-Connector Cable:

- Steam Sterilization
- Gravity Displacement Steam Sterilization
- STERRAD® Sterilization

Steam Sterilization

Prevacuum: Wrapped: 132°C–135°C (270°F–275°F) for 3–4 min.
Unwrapped: “flash” 132°C for 4 min.

Gravity Displacement Steam Sterilization

Wrapped: 132°C–135°C (270°F–275°F) for 15 minutes
Unwrapped: “Flash” 132°C–135°C for 15 minutes

STERRAD® Sterilization

The TRANSDISCAL® Y-Connector Cable may be sterilized with the following STERRAD® systems:

- STERRAD® 100S
- STERRAD 50
- STERRAD 200
- STERRAD NX®
- STERRAD 100NX

All instructions given in the corresponding STERRAD® Sterilization System User's Guide must be followed.

Note: The TRANSDISCAL® Y-Connector Cable should NOT be sterilized within the autoclave case. Any validated tray recommended for use with STERRAD may be used.

⚠ Warning:

Halyard Health has validated ONLY these cleaning and sterilization methods for the TRANSDISCAL® Y-Connector Cable. No other cleaning and sterilization methods have been tested. If any other type of cleaning or sterilization method is used on this product, it is the user's responsibility to verify sterility. Failure to properly clean the device can lead to patient injury.

Troubleshooting

The following table is provided to assist the user in diagnosing potential problems.

PROBLEM	COMMENTS	TROUBLESHOOTING
No temperature measurement in Treatment Mode OR Inaccurate, erratic or sluggish temperature reading in Treatment Mode	In order to measure temperature the entire system must be connected and all devices must be in good working order.	Ensure that all equipment required is correctly connected. Check for an error message on the RF Generator. Visually inspect the probe(s) and cable(s) for damage. Ensure that devices are dry and at room temperature. If problem persists, discontinue use.
No impedance measurement in Placement Mode	In order to measure impedance, the entire system must be connected, and all devices must be in good working order.	Ensure that all equipment required is correctly connected. Check for an error message on the RF Generator. Visually inspect the probe(s) and cable(s) for damage. Ensure that devices are dry and at room temperature. If problem persists, discontinue use.
Probe Connector does not fit in Probe Plug-In	Each connector is designed to connect in a specific way for safety reasons. If the connector “keys” are out of line, the connectors won't fit together.	Check that the connector keys are lined up in the proper orientation. Check that the connectors being connected have the same number of pins and that the pins are not damaged. Ensure that the connectors are clean and unobstructed.

Customer Service and Product Return Information

If you have any problems with or questions about this HALYARD® Equipment, Halyard Health
5405 Windward Parkway
Alpharetta, GA 30004 USA
E-mail: PMPorders@hyh.com
1-844-425-9273 (1-844-HALYARD)

Notes

In order to return products under limited warranty you must have a return authorization number before shipping the products back to Halyard Health.

Limited Warranty

Halyard Health warrants that these products are free from defects in original workmanship and materials. If these products prove to be defective in original workmanship or original materials, Halyard Health, in its absolute and sole discretion, will replace or repair any such product, less charges for transportation and labor costs incidental to inspection, removal or restocking of product.

This limited warranty applies only to original factory delivered products that have been used for their normal and intended uses. Halyard Health's limited warranty shall NOT apply to Halyard Health's products which have been repaired, altered or modified in any way and shall NOT apply to Halyard Health's products which have been improperly stored or improperly installed, operated or maintained contrary to Halyard Health's Instructions. The warranty period for HALYARD® RF Probes and RF Connector Cables is 90 days from the date of purchase, unless otherwise stated.

Disclaimer and Exclusion of Other Warranties

There are no warranties of any kind, which extend beyond the description of the warranties above. Halyard Health disclaims and excludes all warranties, whether expressed or implied, of merchantability or fitness for a particular use of purpose.

Limitation of Liability for Damages

In any claim or lawsuit for damages arising from alleged breach of warranty, breach of contract, negligence, product liability or any other legal or equitable theory, the buyer specifically agrees that Halyard Health shall not be liable for damages for loss of profits or claims of buyer's customers for any such damages. Halyard Health's sole liability for damages shall be limited to the cost to buyer of the specified goods sold by Halyard Health to buyer which give rise to the claim for liability.

The buyer's use of this product shall be deemed acceptance of the terms and conditions of these limited warranties, exclusions, disclaimers and limitations of liability for money damages.

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www.halyardhealth.com

Distributed in the USA by Halyard Sales, LLC, Alpharetta, GA 30004; In USA, please call 1-844-425-9273 • halyardhealth.com

 Halyard Health, Inc., 5405 Windward Parkway, Alpharetta, GA 30004 USA

 Halyard Belgium BVBA, Leonardo Da Vincilaan 1, 1930 Zaventem, Belgium

Sponsored in Australia by Halyard Australia Pty Limited; 52 Alfred Street, Milsons Point, NSW 2061

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*Registered Trademark or Trademark of Halyard Health, Inc., or its affiliates. © 2015 HYH. All rights reserved.

KIMBERLY-CLARK is a Registered Trademark of Kimberly-Clark Worldwide, Inc.

STERRAD® is a Registered Trademark of Johnson & Johnson.

2016-01-25

15-H1-368-0-00 / 70169425