

ENJOY MORE FACE TIME

Seeing a child or loved one in the hospital is distressing enough, but seeing tubes and tape on their face can cause even more distress.

With ANCORIS* you can *escape from tape* and keep nasoenteric feeding tubes simply, securely, and neatly in place—without the hassle, mess, and distress of tape.




Feels better, looks better than tape

- Eliminates the need to cover the cheek or nose with tape or other adhesive tube securement devices.
- Minimal skin contact may reduce the risk of irritation.
- Flexible positioning on the face enables the clinician to secure the tube out of the patient's line of sight.
- Compatible with the most commonly used small-bore nasoenteric feeding tubes on the market today.



PUT A SMILE ON THEIR FACES

Give your patient and their families something to smile about.



ANCORIS* features an innovative clamp and extended-wear adhesive proven to secure nasogastric feeding tubes for an average of 3 days.[†]

Without the mess of tape!



Item No.	Compatible Tube French Size	Case Qty
ATAS-6-NSB	6 FR, 6.5 FR	10
ATAS-8-NSB	8 FR	10
ATAS-10-NSB	10 FR	10
ATAS-12-NSB	12 FR	10
ATAS-14-NSB	14 FR	10


For more information please visit:

www.avanos.com.au

1800 101 021 in Australia

0800 525 270 in New Zealand

[†]Results from clinical evaluation of product in use. Data on file.



AVANOS

avanosmedicaldevices.com



* Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.

©2023 AVNS. All rights reserved. COPY-06379

Availability of these products might vary from country or region as a result of specific local regulatory approval or clearance requirements for sale in such country or region. Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.